



Cardiology Stethoscope

TSphonette

Guidance

Rev.1.0

Table Of Contents

| | |
|--|---------|
| ■ Product Outline | _____3 |
| ■ What is Extensible diaphragm? | _____4 |
| ■ Comparison of Frequency Characteristics | _____5 |
| ■ How to Use (How to press the chest piece) | _____6 |
| ■ A paper was published | _____7 |
| ■ Ear Friendly Design | _____8 |
| ■ High grade model and Standard model | _____9 |
| ■ Q&A | _____10 |

Product Outline

“ONLY ONE” stethoscope

produced with **Clinical × Engineering**
(Dr. Takashina) (Dr. Shimizu)

1) The world's first !!

Wide Range & Powerful Sound

By using “Extensible Diaphragm”

- 2) Ear friendly design such as Movable Ear Piece and Movable Binaural
- 3) Chest Piece and Binaural are made of Stainless-steel
- 4) **High Grade Model** and **Standard Model** are carried

TSphonette *Grande* No.182



Black
Burgundy
Dark Navy

High grade model



- CHP & Binaural with mirror finish
- Two In One Tube
- Molded Outer Spring (10 years warranty for broken spring)

TSphonette No.132

Black
Burgundy
Dark Navy




Standard Model




- CHP with mat finish
- Single Tube
- Inner Spring

What is Extensible Diaphragm ?



Elastic diaphragm expands in the radial direction



Wide range and **Powerful sound** are generated by transducing repulsion which generate by stretching tight of this elastic diaphragm into sound energy.

※The conventional diaphragm simply uses the vibration characteristics of resin.

The extension diaphragm consists of an extensible elastic membrane part and a highly rigid transparent film part. When the diaphragm is fixed with a ring, the inner ring provided inside the chest piece causes the diaphragm to extend radially over the entire circumference.

Comparison of Frequency Characteristics




The above chart shows comparison of frequency characteristics of SDS and EDS. You can see that the low frequency sound can be heard more clearly with EDS than with SDS.

It is particularly useful for cardiovascular diagnosis because even heart murmur in super low frequency range, such as grade III or IV, can be heard clearly.

How to use (How to press the chest piece)

When pressing the diaphragm lightly into a chest, Every sounds from in high frequency to low frequency are heard in a balanced manner. When pressing more firmly, low frequency sound will be intensified and you can hear every sounds more strongly and even heart murmur Grade III and IV can be heard clearly.

※By changing the pressing force, the elastic membrane part expands and this make frequency characteristics change.



A paper was published

Japanese Circulation Society official publication

『**Circulation Journal** vol.80 No.9 September2016 (P.2047-P.2049) 』



Outline

“This study compared the diagnostic efficacy of the common suspended diaphragm stethoscope (SDS) with a new extensible diaphragm stethoscope (EDS) for low-frequency heart sounds. Based on the results of the sound analysis, the EDS is more efficient than the SDS.”

Developer;
 Dr. Tsunekazu Takashina Head Director, Japan Educational Clinical Cardiology Society
 Dr. Masashi Shimizu emeritus professor, Tokyo Institute of Technology

Ear Friendly Design

Selectable Ear Piece

You can use your favorite one

Movable Soft Ear Piece



Sealing Ear Piece




Ear Piece Hard type



Movable Binaural

You can adjust to your comfortable angle



Tough spring of moderate pressure

You can wear it without ear pain.


No.182

Molded Outer Spring
Pressure of spring is adjustable.
10 year warranty for broken spring



No.132

Inner Spring
You can wear it with moderate pressure



High Grade Model & Standard Model

※Chest piece and Binaural of both model are made of stainless-steel which have better acoustic property

(Mainly for doctors)

High Grade Model

TSphonette *Grande* No.182

Binaural

Mirror finish with premier feel
Made of stainless-steel

Molded Outer Spring

Pressure of spring is adjustable.
[10 year warranty for broken spring](#)



Chest Piece (CHP)

Mirror finish with premier feel
Made of stainless-steel

Two In One Tube

For better acoustic property



3 colors
• Black
• Burgundy
• Dark navy



(Mainly for nurses)

Standard Model

TSphonette No.132

Binaural

made of stainless-steel

Inner Spring

You can wear it with moderate pressure



Chest Piece (CHP)

Mat finish by special polishing
Made of stainless-steel

Single Tube

3 colors
• Black
• Burgundy
• Dark navy



High Grade Model & Standard Model

High Grade Model

Standard Model

| Product Name | Tsphonette Grande No.182 | TSphonette No.132 |
|------------------------|--|----------------------------|
| User | Mainly doctors | Mainly nurse |
| Length | About 710mm | About 710mm |
| Weight | About 230g | About 160 g |
| Diaphragm | Φ52mm | Φ52mm |
| Material(CHP/Binaural) | Stainless-steel | Stainless-steel |
| Finishing | Mirror finish (CHP/Binaural) | Mat finish (CHP) |
| Binaural, Spring | Φ6mm · Molded outer spring | Φ5mm · Inner spring |
| Tube | Two in one tube | Single tube |
| Color | Black, Burgundy, Dark Navy | Black, Burgundy, Dark Navy |
| Warranty | One year (10 years warranty for broken spring) | One year |
| | | |